


## MaxTwo Neuronal Cell Plating Protocol

### 1. Experiment Workflow



\* These are recommended, however other experiments may require specific frequency of media changes and timing of neural activity recording.

### Notes

Change 50% medium twice a week from Day 1\*.

### 2. Required Materials

| Item                                      | Supplier                          | Catalog Number                 |
|---|-----------------------------------|--------------------------------|
| MaxTwo Multi-Well Plate                   | MaxWell Biosystems AG             | MX2-S-6W                       |
| Frozen Primary Rat Neurons                | Brain Bits*                       | FSDECX1M                       |
| Human iPSC-derived Neurons                | BrainXell*<br>Elixirgen*<br>FCDI* | BX-0500<br>EXGS-QNMSV<br>R1088 |
| Cell Culture Medium                       | Multiple Vendors                  |                                |
| Ethanol 70%                               | Multiple Vendors                  |                                |
| Sterile Deionized Water                   | Multiple Vendors                  |                                |
| Surface Coating Materials (see section 5) | Multiple Vendors                  |                                |
| Terg-a-zyme                               | Sigma-Aldrich                     | Z273287                        |

\* These are example sources of neuronal cells that we have tested in MaxOne and MaxTwo. If you are interested to establish a protocol for your cells, contact us at [info@mxwbio.com](mailto:info@mxwbio.com).

### 3. MaxTwo Multi-Well Plate Sterilization

#### Procedure

1. Prepare Terg-a-zyme-solution in deionized water (1% solution: 10 g per liter).
2. Fill each well of the MaxTwo Multi-Well Plate with 1 mL of 1% Terg-a-zyme-solution and leave over night at room temperature (R.T).
3. Wash wells three times with deionized water.
4. Remove the remaining sterile deionized water and spray the MaxTwo Multi-Well Plate and MaxTwo lid thoroughly with 70% ethanol (see Figure 1).
5. Fill each well and each compartment of the MaxTwo Multi-Well Plate with 70% ethanol (see Figure 2).
6. Transfer the MaxTwo Multi-Well Plate to the biological safety cabinet.
7. Remove the 70% ethanol from the wells and compartments after 30 minutes.
8. Rinse each well and compartment of the MaxTwo Multi-Well Plate 3 times with sterile deionized water.
9. Aspirate the sterile water with a vacuum pump. **[A]**
10. Fill the small compartments between the wells with 0.5 mL and the big compartments with 1 mL sterile deionized water (see Figure 3).
11. Cover the MaxTwo Multi-Well Plate with its Lid.
12. Proceed to surface coating.

**A:** Make sure that the bottom of the MaxTwo Multi-Well Plate is dry!




**Figure 1:** Sterilize the MaxTwo Multi-Well Plate



**Figure 2:** Fill the wells of MaxTwo Multi-Well Plate with 70% ethanol



**Figure 3:** Fill compartments between wells with sterile deionized water



## 4. MaxTwo Multi-Well Plate Surface Coating

Notes

### Primary Coating Chemical Options

|                       |   |   |                                      |
|-----------------------|---|---|--------------------------------------|
| Chemical              | Poly (Ethyleneimine) (PEI)  | Poly-L-Ornithine (PLO)  | Poly-D-Lysine (PDL)                  |
| Company               | Sigma/Merck   | Sigma/Merck   | ThermoFisher Scientific              |
| Catalog Number        | P3143   | P4957   | A3890401                             |
| Working Concentration | 0.07%   | 0.005%  | 0.1 mg/mL                            |
| Incubation Location   | 37°C<br>5% CO <sub>2</sub> incubator                                  | 37°C<br>5% CO <sub>2</sub> incubator                                    | 37°C<br>5% CO <sub>2</sub> incubator |
| Incubation Time       | 1 hour  | 2 hours   | 1 hour                               |
| Preparation           | Dilute stock in sterile borate buffer<br>See preparation notes<br>PEI | Dilute stock in sterile deionized water<br>See preparation notes<br>PLO | Ready to use                         |

#### Preparation for PEI


1. Prepare 100 mL of 1X borate buffer: dilute 5 mL of 20X borate buffer (ThermoFisher Scientific 28341) in 95 mL in sterile water.
2. Prepare an intermediate 7% PEI solution: pour 1 mL of 50% PEI solution into a 15 mL centrifuge tube and allow it to settle. Add 6 mL of 1X borate buffer to obtain an intermediate ~7% PEI solution (can be stored in 1 mL aliquots at -20°C for 1 month).
3. Prepare a final ~0.07% PEI solution by diluting 1 mL of intermediate 7% PEI solution in 99 mL of 1X borate buffer.
4. For sterilization, filter through a 0.22 µm filter unit.

#### Preparation for PLO

1. Mix 1 mL 0.01% PLO stock solution with 1 mL sterile water.
2. For sterilization, filter through a 0.22 µm filter unit.

## Surface Coating Area Options

## Notes

|   | Whole-Area Plating  |
|---|---|
| Area                                      | 25 mm <sup>2</sup>  |
| Drop for primary coating (e.g. PDL)       | 50 µL   |
| Drop for secondary coating (e.g. Geltrex) | 50 µL   |
| Cell drop volume                          | 50 µL   |
| Primary rat neurons                       |   |
| Cell number                               | 60'000 - 80'000   |
| Density (cells/mm <sup>2</sup> )          | 2'400 - 3'200   |
| Human iPSC-derived neurons                |   |
| Cell number                               | 200'000 - 300'000   |
| Density (cells/mm <sup>2</sup> )          | 8'000 - 12'000  |
| Picture of plating area                   |  |

### Procedure (example PDL)

1. Add PDL solution to each well of the MaxTwo Multi-Well Plate covering the entire electrode array (see Figure 4). **[B]**
2. Cover the MaxTwo Multi-Well Plate with its Lid.
3. Incubate the MaxTwo Multi-Well Plate with PDL in a 5% CO<sub>2</sub> incubator at 37°C. **[C]**
4. Aspirate PDL.
5. Rinse each well of the MaxTwo Multi-Well Plate 3 times with 1 mL sterile deionized water.
6. Aspirate the sterile water with a vacuum pump.

**B:** Gently pipette the coating solution onto the center of each MEA. Avoid touching the MEA surface with the tip of the pipette!

**C:** Adapt incubation time according to primary coating chemical.

## 5. Cell Plating

### Secondary Coating Chemical Options


| Chemical              | Laminin                              | Matrigel                             | Geltrex                              |
|-----------------------|--------------------------------------|--------------------------------------|--------------------------------------|
| Company               | Sigma/Merck                          | Corning                              | ThermoFisher Scientific              |
| Catalog Number        | L2020                                | 354230                               | A1569601                             |
| Working Concentration | 0.02 mg/mL                           | 0.4 mg/mL                            | 1x                                   |
| Incubation Location   | 37°C<br>5% CO <sub>2</sub> incubator | 37°C<br>5% CO <sub>2</sub> incubator | 37°C<br>5% CO <sub>2</sub> incubator |
| Incubation Time       | 30 minutes                           | 2 hours                              | 1 hour                               |
| Preparation           | Dilute stock in cell culture media   | Dilute stock in cell culture media   | Ready to use                         |

### Procedure (example Geltrex)

1. Add Geltrex solution to the center of each well of the MaxTwo Multi-Well Plate (see Figure 4).
2. Cover the MaxTwo Multi-Well Plate with its Lid.
3. Incubate the MaxTwo Multi-Well Plate with Geltrex in a 5% CO<sub>2</sub> incubator at 37°C. **[D]**
4. Aspirate Geltrex.
5. Add cell solution to each well of the MaxTwo Multi-Well Plate on the area previously covered by Geltrex. **[E]**
6. Cover the MaxTwo Multi-Well Plate with its Lid.
7. Incubate MaxTwo Multi-Well Plate with neurons in a 5% CO<sub>2</sub> incubator at 37°C for 1 hour. **[F]**
8. After incubation, fill up each MaxTwo Multi-Well Plate carefully with 1.2 mL of medium.
9. Change 50% of the medium on day 1 after plating.
10. Maintain the cell cultures inside a 5% CO<sub>2</sub> incubator at 37°C to control environmental conditions. Replace 50% of the medium twice a week.

### Notes


**D:** Adapt incubation time according to secondary coating chemical.

**E:** Make sure to mix regularly the tube with cells when dotting several wells, otherwise cells will settle at the bottom of the tube. 

As much as possible, keep the cell droplet at the center of the electrode array.

**F:** Covering the MaxTwo Multi-Well Plate is crucial to avoid medium evaporation while cells settle down on the surface.

**Figure 4:** Placing coating solution



#### Conditions of Use

For life science research use only.

#### Revision History

Version 1: November 2020

© MaxWell Biosystems AG. All rights reserved.