

MaxTwo

Multiwell HD-MEA System

Maximize Your Cell Functional Assays



Every Cell Has a Story to Tell


Next-generation multiwell high-density microelectrode array system for functional imaging of your cells.

maxwell
BIOSYSTEMS


High-Density Microelectrode Array (HD-MEA)

MaxTwo


Capturing the intricate dynamics of cellular activity is key for various applications, including phenotype characterization, disease modeling, and drug discovery. With the rapid advancements in in-vitro cell and tissue models, there is an increasing demand for real-time and label-free methods to simultaneously study the function, maturation, connectivity and morphology of cells within the same sample. High-density microelectrode arrays (HD-MEAs) serve as an unparalleled tool for conducting non-invasive in-vitro functional assays.



MaxTwo HD-MEAs represent the next-generation electrophysiology platform, featuring densely packed microelectrodes capable of capturing the electrical signature of electrogenic cells with unprecedented precision.



26'400 Pt-Electrodes Per Well
3'265 Electrodes/mm²
Low-Noise Readouts
Flexible Electrical Stimulation



Readout at different scales
Network level
Cell level
Subcellular level


Featuring 26'400 electrodes within a large sensor area of 2x4 mm² per well, MaxTwo HD-MEAs measure the activity of cells across multiple scales – capturing whole network dynamics and identifying single cell function. MaxTwo HD-MEAs' high electrode resolution and low noise ensures that spikes generated by any neuron on the array can be picked up by an electrode in close vicinity. This excellent sensitivity enables to detect cells' small signals, such as in developing induced pluripotent stem cell (iPSC)-derived neurons and propagating action potentials along axonal arbors. MaxTwo HD-MEAs' exceptional data quality enables you to derive fast and meaningful conclusions in your cell-based assays.

Key Features




High-Quality Data

Obtain high-resolution and high-quality electrophysiological data while tracking dynamic functional changes at cellular, subcellular, and network levels.




Sensitivity

Capture the smallest signals (μV) thanks to our low-noise recording channels and the high electrode density, ensuring an electrode in close vicinity to even the weakest signals.



Throughput

Record from 6 or 24 wells, increasing testing capabilities and decreasing experimental time, while ensuring consistency across wells.



Reproducibility

Exploit optimized recording strategies to analyze the entire culture at individual neuronal levels, thus enhancing the data reproducibility and statistical power.



Long-Term Experiments

Assess cell development, maturation, or compound effects by performing longitudinal experiments over the course of days, weeks, and months.



Integrity

Conduct non-invasive and label-free recordings, eliminating any potential side effects associated with the use of dyes and prolonged exposure to light.

Applications

Exploring the electrophysiological properties of cells enables to understand functions in healthy and disease states. Additionally, positive or adverse compound effects can be revealed by changes in electrical signatures of cells and networks. An extracellular electrophysiology platform like MaxTwo allows to investigate disease progression, phenotype variations, and the impact of compounds over extended time, all within a multiwell format.



Phenotyping

Solution for extended cell characterization and longitudinal developmental studies spanning weeks or months.

Disease Modeling

Versatile platform to study various cellular models of neurodevelopmental and neurodegenerative diseases.

Drug Discovery

Monitoring of cellular responses to compound administration, ensuring a comprehensive evaluation of drug efficacy.

Toxicology

Ideal tool to assess safety and toxicity levels of compounds in different cellular models.

Assays


MaxLab Live Software is a comprehensive and user-friendly tool designed for visualizing, recording, and analyzing data recorded with MaxTwo. Its assays guide the user through automated workflows to the extraction of key metrics at multiple levels, including network, single-cell, and subcellular details.

ActivityScan

Are the cells active?
Where are they located?

Record whole sample activity through array scanning

Identify the firing activity and location of cells



Activity Metrics


Firing Rate, Spike Amplitude, etc.

Network

Are the cells interconnected?

Smart electrode selection to record population activity

Identify network-level synchronized activity



Network Metrics


Burst Rate, Burst Duration, etc.

AxonTracking

Are the axons functional?

Record and automatically identify single neurons

Track axonal arbors and extract action potential conduction velocity



Axon Metrics


Conduction Velocity, Total Axonal Length, etc.

Stimulation

How does electrical stimulation affect cells?

Characterize the excitability of cells


Study stimulation effects on neuronal plasticity



Stimulation Metrics



Evoked Total Spikes, Evoked Activity Peak, etc.

Experimental Workflow




Case Study

24-Well Plate



ALS Disease Modeling with iCell Motor Neurons: Control (AHN) and disease models TDP43 (M337V), TDP43 (Q331K). A), B), D) show data recorded at days in vitro (DIV) 21 from cells co-cultured with iCell Astrocytes 2.0 (C1249) on a MaxTwo 24-Well Plate. A) Raster plot of each well highlights well-to-well consistency (conditions grouped by column). B) Representative Network Activity plots (left to right: AHN, M337V and Q331K). Please note the y axis difference. C) Burst duration plotted over 4 weeks for AHN, M337V and Q331K. D) Example Axon Traces (left, AHN; right, M337V). E) Total Axon Length over 4 weeks for AHN, M337V and Q331K. Data courtesy: Fujifilm Cellular Dynamics, Inc.

MaxTwo System



MaxTwo System Features

One-button interface	Open/close the cover Turn on/off the system Lock the well plate
System status indicator	LED
Dimension (LxWxH)	40x16x12 cm ³
Well-Plate compatibility	6- and 24-Well Plate
Temperature control	Yes
CO ₂ control	Yes, via gas mixer

MaxTwo 6-Well Plate



MaxTwo 24-Well Plate



MaxTwo Plate Features

6-Well

24-Well

MaxTwo Plate Features	6-Well	24-Well
Well-Plate dimension (L x W x H)	127.8 x 85.8 x 16.1 mm ³	127.8 x 85.8 x 16.1 mm ³
Number of wells	6	24
Well (inner) diameter	29.0 mm	15.3 mm
Active sensing area	3.85 x 2.10 mm ²	3.85 x 2.10 mm ²
Total number of electrodes	6 x 26'400	24 x 26'400
Electrode density	3'265 electrodes/mm ²	3'265 electrodes/mm ²
Electrode center-to-center distance	17.5 μm	17.5 μm
Electrode size	12.0 x 8.8 μm ²	12.5 x 12.5 μm ²
Electrode material	Platinum black	Platinum
Surface topography	<2.0 μm	<0.3 μm
Sampling rate	10.0 kHz/channel	10.0 kHz/channel
Number of recording channels	1'020/well	1'020/well
Electrical stimulation	Yes	No
Simultaneous recordings	All 6 wells in parallel	6 wells in parallel, row by row
Sterilization	No	EtO treatment according to EN ISO 11737-2

Disclaimer

For life science research use only. The contents of this document are provided by MaxWell Biosystems, 'as is'. MaxWell Biosystems makes no representations nor warranties with respect to the accuracy or completeness of the contents of this publication and reserves the right to make changes to the specification at any time without notice. All trademarks are the property of their respective owners.

maxwell

BIOSYSTEMS


 www.mxwbio.com

 MaxWell Biosystems

 @mxwbio

 info@mxwbio.com

 +41 44 244 24 24

 Albisriederstrasse 253
8047 Zurich, Switzerland

